As men age, millions suffer from the symptoms of BPH (benign prostatic hyperplasia), also known as enlarged prostate. The good news is that BPH is not life threatening and can be treated. Some of the symptoms of BPH are:

- Frequent and urgent need to urinate
- Weak flow
- Inability to completely empty the bladder
- Sleep disruption as a result of needing to urinate frequently at night

While some men may feel embarrassed to talk about it, BPH is a significant quality of life issue that deserves frank discussion. The U.S. Food and Drug Administration has approved a device called the Prolieve™ Thermodilatation System for the treatment of BPH symptoms. Treatment with this device:

- is safe and effective
- provides significant symptomatic relief
- is a minimally invasive procedure that can be performed in a single 45 minute office visit
- offers a viable alternative to treatment by medications or surgery
- 99% of patients did NOT experience any form of erectile dysfunction following treatment
- 94% of patients reported no to mild pain during treatment
- 82% of patients go home without a catheter

While BPH typically does not pose a major health risk, it can affect a man's quality of life and often significantly impedes even the simplest daily activities. It is typical for men and their partners to feel tired due to the sleep disruption caused by the man's need to urinate frequently at night. At the same time some men may not even realize they have a treatable condition.

How does the Prolieve system work?

The physician will insert a catheter into your urethra. When microwave energy is turned on, it will be transmitted from a special antenna inside the catheter to provide continuous heat to your prostate. A small balloon that is also part of the catheter will inflate within the section of the urethra close to where the prostate is located. The inflated balloon will cause the tissue to expand. During the entire 45 minute procedure, warm water will be circulating inside the catheter through your urethra. For your safety, the Prolieve System's computer monitors the temperature surrounding the treatment area by means of a rectal temperature monitor. If the temperature reaches 42 C, the system will shut off automatically.

The Prolieve System is FDA approved. The System's design and operating protocol were investigated in clinical studies of other men suffering from BPH. The results showed the Prolieve System is safe and effective.

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